



Portable Dental X Ray(**YJ-DXP02**) and sensor(**YJ-HDR500A**) solution

2. Product description:

Features(**YJ-DXP02**):

- 1.Computer Touch control.
- 2.Toshiba D-045 Tube.
- 3.High Frequency DC.
- 4.Replaceable battery.



Parameters:

Model	YJ-DXP02
X Ray Tube	Toshiba D-045
Output Voltage/Current	65kv/2mA
Focus Spot	0.4mm
Anode Angle	16



Forever Medical

Dental x-ray More picture:



Features of Dental X Ray Sensor YJ-HDR500A

Long lifespan, high definition, high efficiency, high benefit, Dicom.

Multiple user solution. Powerful user solution. Powerful image process software.



Specifications

Model	YJ-HDR500A
Chip Type	CMO APS
Fiber plates	YES
Scintillator	GOS
Size:	Size 1 :27.5*38.5
Effective area	Size 1 :22.5*30
Pixel dimension	Size 1 :375mm3
Pisel	18.5μm
Resolution	Size 1:1.9m(1600*1200)
Operating system	Theore tical:27
Others	True:>20
Multi-language	Czech,English,French,German,Hungarian,Italian,Polish,Portuguese,Russian,Spanis h,Turkish,Chinese.



Forever Medical

Dental sensor image :

